www.ocivm.com

Vacuum Microengineering Inc. jociepa@ocivm.com

## Thermal diffusion behaviour of Li on surfaces of HOPG, Si(111) and SiC-6H measured with LEED and AES

Paige Harford, Jozef Ociepa, Guenevere O'Hara, Bart Checinski OCI Vacuum Microengineering, Inc., 200 Stronach Crescent, London, ON Canada

## INTRODUCTION

- **GOAL:** To observe lithium reactivity and mobility on highly oriented pyrolytic graphite (HOPG), silicon (111), and silicon carbide (6H)
  - Samples are potential anode and cathode materials for lithium-ion batteries
- Use of Low Energy Electron Diffraction (LEED) and Auger Electron Spectroscopy (AES) for analysis

## MOTIVATION


Most current battery development research is focused on the improvement of battery properties such as:

- Maximum capacity
- Power density
- Solid state electrolyte solutions
- Cycle life

In this study, we direct our attention to fundamental aspects of the battery anode materials with respect to lithium diffusion. Deposition of Li on surfaces of relevant materials with structural characterization will demonstrate thermal Li behaviour, presenting possible links to Li diffusion mechanisms for novel battery designs.


## **EXPERIMENTAL SETUP**

Testing was completed in an Ultra-High Vacuum (UHV) system utilizing Li evaporation and high-temperature annealing (Model IMBE300-SPH). Li evaporation was monitored with a quartz crystal microbalance.



# Significant room temperature Li diffusion on HOPG and SiC-6H associated with unchanged surface crystallography

Before Lithium Evaporation



## No room temperature Li diffusion on Si(111) observed in drastic surface crystallography change

ISCUSSION

The reactivity of Si makes it a challenging anode material, despite high theoretical capacity


Results depicted a persistent AES Li peak during annealing and crystallography destruction immediately following evaporation (evident in LEED data)

Literature shows interstitial Li diffusion with interactions between native Si atoms and Li atoms. Interaction causes expansion and destroys surface crystallography (as shown in LEED data)




Future battery development could benefit from this methodology of using AES to observe thermal lithium diffusion tendencies for characterization of possible solid-state battery materials.


5Å Lithium Evaporation



The LEED surface crystallography data and AES data taken from the sample cover and sample plate on both samples confirms the diffusion of lithium at room temperature, without impairing surface crystallography

Literature shows interstitial Li diffusion as a mechanism in HOPG and SiC-6H; however, LEED and AES data suggests native atom interaction is not crystallography-altering






Model MBE300-GB



### 3 Annealed from 50°C to 800°C



→ Room Temperature Lithium Deposition

- AES scans on plate confirm successful Li deposition on HOPG and SiC-6H
- Timed AES scans on SiC-6H reveal room temperature diffusion over the span of approximately 1 hour

## **FUTURE RESEARCH**

- Probing anode and cathode materials with a liquid electrolyte
  - Following Li deposition and diffusion, the sample can be moved to the liquid electrolyte (in the glovebox) and moved back after interactions for further testing (LEED/AES characterization before and after electrolyte contact)
  - Full battery integration is also possible in this configuration